medel











Veloce / Fast
Facile / Easy

NO CONTACT

TERMOMETRO INFRAROSSI A DISTANZA / TOUCH FREE INFRARED THERMOMETER

IT

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INTRODUCTION

Utilizing infrared technology, this thermometer takes temperatures in seconds by measuring heat generated by the surface skin of the forehead or other objects.
 This product conforms to the provisions of the EC directive MDD/93/42/EEC.

Its advantages include:

1. 6 in 1 Functions: Human Body/Object/Night Mode Option/Fever Alarm/30 User Memory/Backlight/2 Night Mode Option: Switch to Night mode to reduce the interference of buzzer during your baby sleeping.
3. Red LED Light for Fever Alarm: This device has the LED light to remind user their readings are over 38°C and plus the beep sound only in Human Body Mode.
4. 30 User Memory
5. Illuminated Backlight Display
6. °C/°F Switch-able Function
7. One-second Reading
8. Auto power off for power saving
9. Low-battery Indicator: Indications for battery condition and measuring range.
10. Large LCD Display

IMPORTANT INFORMATION BEFORE USE

When using this product, please be sure to follow all the notes listed below. Any action against these notices may cause injury or affect the accuracy.

1. Do not disassemble, repair, or remodel the thermometer.
2. Be sure to clean the thermometer lens each time after usage.
3. Avoid direct finger contact with the lens.
4. No modification of this equipment is allowed.
5. It is recommended that user may take 3 temperatures. If they are different, use the highest reading.
6. Do not expose the thermometer to extreme temperature, very high humidity, or direct sunlight.
7. Avoid extreme shock or dropping the device.
8. Before the measurement, patients and thermometer should stay in steady state room condition for at least 30 minutes.
9. Avoid measuring temperature in 30 minutes after exercise, bathing, or returning from outdoor.
10. To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.
11. It is ill-advised to disassemble the thermometer.
12. Please use the thermometer solely for its intended purpose.
13. Carefully hold the device when in use to avoid dropping the device.
14. Allow one minute between successive measurements as slight variations may occur if measurements are taken over a short period of time. Use average temperatures instead.
15. There is no absolute body temperature standards. Keep reliable records of your personal temperature to serve as a reference for judging a fever.
16. Under any circumstances, the temperature taking result is ONLY for reference. Before taking any medical action, please consult your physician.

PRODUCT IDENTIFICATION (FIG. A)

1. LCD Display
2. Red LED Light for Fever Alarm (Only in Human Mode)
3. Human/Object Button
4. Power/Measure button
5. Memory Button
6. Measurement Sensor
7. Buzzer
8. Battery Compartment lid

DESCRIPTION OF LCD DISPLAY (FIG. B)

1. MEASURE in Progress
2. Night Mode
3. Human MEASURE Mode
4. Object MEASURE Mode
5. Degrees of Temperature
6. Celsius Scale
7. Fahrenheit Scale
8. Low battery Warning
9. Smiling Face (Indicating temperature is lower than 38 °C or 100.4°F)
10. Frowning Face (Indicating temperature is equal or higher than 38°C or 100.4°F)
11. Memory Symbol and Sets of Records

BATTERY INSTALLATION

Low battery warning:
 When the battery power becomes low, the low battery symbol "Lo" will appear on the display. The thermometer can still be used during this time, but the batteries should be replaced as soon as possible. If the batteries run out completely, "Lo" will be displayed along with the low battery symbol. In this case,

the batteries will need to be replaced before using the thermometer again. (FIG. C)

Replacing the Battery:

1. Gently slide the battery cover back.
 2. Carefully remove the old batteries and properly discard.
 3. Insert new batteries (Two 1.5V alkaline AAA Size) according to the proper polarity. (FIG. D)
 4. Slide the battery cover back on.
- NOTE: Battery-operated**
1. Please properly dispose of the batteries away from small children and heat.
 2. It is recommended to remove the batteries if the unit will not be used for an extended period of time.
 3. For long durations of non-operation, please remove all batteries from the device.
 4. Batteries must be disposed of in accordance with local environmental and institutional policies.
 5. Dispose of used batteries in accordance with the applicable legal regulations. Never dispose of batteries in the normal household waste

SWITCHING BETWEEN FAHRENHEIT OR CELSIUS

Your thermometer can display results in either degrees Celsius (°C) or degrees Fahrenheit (°F). To switch between Celsius and Fahrenheit, while the unit is on press and hold both the Human/Object button and the Memory button for approximately 3 seconds. This will change the mode to either °C or °F. (FIG. E)

Switching between 4 Kinds of MEASURE Mode

1. Under power on status, you can press the Human/Object button (MODE) to switch different MEASURE mode. There are 4 kinds of mode which including Human, Object, Human/Night, and Object/Night mode (in order).
 2. The beep sounds will be closed when your choice in Human/Night mode, and Object/Night mode, and the Moon symbol will appear on the LCD in both Night mode. (FIG. F)
 a. Human Mode b. Object Mode
 c. Human Mode & Night Mode d. Object Mode & Night Mode
 [NOTE: Each press will come with a beep sound to ensure the setting is activated. (Except both Night modes)]

TIPS FOR MEASURING HUMAN TEMPERATURE

Bear in mind that the thermometer needs to have been in the room in which the measurement is taken for at least 30 minutes before use.

- 1. Attempting to take temperature readings from sites on the body other than the forehead may produce inaccurate results.
- 2. The patient should remain still while the reading is being taken.
- 3. Infrared forehead temperature readings are equivalent to oral temperature readings. In all of these cases, please consult your doctor.
- 4. Readings taken while asleep should not be compared directly to readings taken while awake, as body temperature while asleep is typically lower.
- 5. Do not take body temperature readings within 30 minutes of being outdoors, exercising or bathing.

MEASURING HUMAN TEMPERATURE (FIG. G)

1. Press the Power button (MEASURE) to turn the thermometer on. The unit will run a self-test and all symbols on the display will momentarily appear.
2. Ensure that the thermometer is in Human mode; the Human symbol will be on the display. To alternate between modes press and release the Human/Object Mode button until you see the desired measurement symbol on the display.
3. Press and hold the "Measure button" (MEASURE), aiming the Measurement Sensor at the person's forehead, holding the thermometer approximately 0.5-3 cm. But don't touch the forehead.
4. Release the "Measure button". In about a second, you will hear a short beep means this temperature reading has been completed and accompany with a back-light.
5. If the temperature measurement is below 38°C, a "Smiling Face ☺" will appear next to the reading. If the reading is 38°C or above, a "Frowning Face ☹" will be displayed and the RED LED light up.
6. After about 30 seconds after use, the thermometer will automatically beep and shut off.

MEASURING OBJECT/LIQUID TEMPERATURE

1. Press the Power button (MEASURE) to turn the thermometer on. All symbols on the display will momentarily appear.
2. Ensure that the thermometer is in Object mode; the Object symbol will be on the display "☾". To alternate between modes press and release the Human/Object Mode button until you see the desired measurement symbol on the display.
3. Press and hold the "Measure button" (MEASURE), aiming the Measurement Sensor at the object, holding the thermometer approximately 0.5-3 cm.
4. Release the "Measure button" (MEASURE) and the temperature reading will be displayed.
5. After about 30 seconds after use, the thermometer will automatically beep and shut off.

MEMORY FUNCTION

Memory Recall:
 You can recall up to 30 measurements currently stored in memory to share with your physician or trained healthcare professional.
 1. When the device is on, press once briefly on the "MEM" button, then pass it again to show the last measurement accompanied by "M" symbol. (Fig. I)
 2. The "M" symbol or "sM" symbol will appear with each measurement stored in memory to indicate whether a person or object temperature was taken.
 3. Each press of the same button recalls a previous measurement, so "M" then all the way to "M" (Fig. J)
Memory Deletion:
 1. Under power on status, you may keep press MEM Button for more than 3 seconds to delete all the readings.
 2. You will see "-" on display and "sound four short beeps", that mean all the memories are cleared. (Fig. K)
 3. Automatically on the 31th measurement; when the 30 memories have been used up, any new measurement will be recorded with "M" and the oldest memory deleted without you having to do anything.
 [NOTE: All the readings will be cleared no matter record in Human mode or Object mode.]

CARE AND MAINTENANCE

Lens/ Measurement Sensor: Gently clean with an alcohol swab. Do not use water to wash the thermometer lens directly.
Thermometer: Clean with a soft, dry cloth. Do not use water to rinse the device

APPLIED STANDARDS

This product conforms to the provisions of the EC directive MDD (93/42/EEC). The following standards apply to design and/or manufacture of the products:
EN 60601-1-2 electromagnetic compatibility;
EN 60601-2-56 Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement
EN 60601-1-2 electric safety standard:
 • Internally powered equipment
 • IP22 (protection against solid foreign objects and against harmful effects due to the ingress of water).
 • Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide
 • Continuous operation

ERROR CODES

LCD Display	Cause	Solution
H i	The temperature measured is higher than: 1. Human thermometer mode: 43°C (109.4°F) 2. Object thermometer mode: 100°C (212.0°F)	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Services.
Lo	The temperature measured is lower than: 1. Human thermometer mode: 34°C (93.2°F) 2. Object thermometer mode: 0°C (32.0°F)	
Err	The operating temperature is not in the range from 15°C-35°C (59°F-95°F)	Operate the thermometer only between the specified temperature ranges.

MAINTENANCE AND REPAIRS

In case of failure, address to qualified personnel authorized by Medel International Srl. Do not open the device in any case. The unit has no user-serviceable parts within and does not need internal maintenance or lubrication.

TECHNICAL SPECIFICATION

MODEL: MEDEL NO CONTACT - TOUCH FREE INFRARED THERMOMETER MODEL NT17
Measuring range : Human Body: 34°C-43°C (93.2°F-109.4°F) Object: 0°C-100°C (32.0°F-212.0°F)
Calibration Accuracy : 35°C-42°C (95°F-107.6°F) ± 0.2°C else : ± 0.3°C
Object : 20°C ± 1°C > 20°C : ± 5%
Display resolution : 0.1°C
Operating environment: 15°C-35°C (59°F-95°F) with relative humidity up to 85% (non condensing)
Storage/ Transportation environment : -25°C-+55°C (-13°F-+131°F) with relative humidity up to 95% (non condensing)
Power supply : 2 x 1.5V AAA size alkaline batteries
Weight : approx. 80g (with batteries)
Dimensions : approx. 128.5mmx48.83mmx38.85mm (LxWxH)
Operation Distance : 0.5-3 cm

EMC TABLES

Guidance and manufacturer's declaration -electromagnetic emissions		
NO CONTACT Thermometer (NT17) is intended for use in the electromagnetic environment specified below. The customer or the user of NO CONTACT Thermometer (NT17) should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment -guidance
RF emissions CISPR 11	Group 1	Medel NO CONTACT Thermometer (NT17) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Medel NO CONTACT Thermometer (NT17) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration -electromagnetic immunity

Guidance and manufacturer's declaration -electromagnetic immunity			
Medel NO CONTACT Thermometer (NT17) is intended for use in the electromagnetic environment specified below. The customer or the user of NO CONTACT Thermometer (NT17) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Recommended separation distances between portable and mobile RF communications equipment and Medel NO CONTACT Thermometer (NT17)

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1.2 vP	80 MHz to 800 MHz d=1.2 vP	800 MHz to 2.5 GHz d=2.3 vP
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
 NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration -electromagnetic immunity

Guidance and manufacturer's declaration -electromagnetic immunity			
Medel NO CONTACT Thermometer (NT17) is intended for use in the electromagnetic environment specified below. The customer or the user of Medel NO CONTACT Thermometer (NT17) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Radiated RF IEC 61000-4-3	80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of Medel NO CONTACT Thermometer (NT17) including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d=1.2 vP 80 MHz to 800 MHz d=2.3 vP 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: ⚡

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Medel NO CONTACT Thermometer (NT17) is used exceeds the applicable RF compliance level above, Medel NO CONTACT Thermometer (NT17) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Medel NO CONTACT Thermometer (NT17).
 b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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WARRANTY CONDITIONS

The device is covered by a 2 years warranty from the date of purchase for any material or manufacturing original defect. The warranty consists in the replacement and/or free repair of originally defective components. The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear. The device must be repaired by authorized technical service centres only. Transport costs for the device are at user's charge. Any repair out of warranty conditions is at user's charge. The warranty is void if the device has been tampered, if the defect results from improper use or if the damage is not imputable to the manufacturer. The warranty does not involve any compensation for any kind of damages, either direct or indirect, to persons or things, occurred when the device is malfunctioning. The warranty is valid from the date of purchase certified by the purchase ticket or invoice which must be peremptorily attached to the warranty coupon.

The lack of the warranty coupon properly filled in and validated by a purchase certification voids the warranty.

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COUPON TO BE RETURNED IN CASE OF REPAIRS THE WARRANTY IS VALID ONLY IF ACCOMPANIED BY THE PURCHASE TICKET/INVOICE

PRODUCT TYPE
 Model: _____
 Serie No.: _____
 Date of purchase: _____
Buyer's DATA
 Name: _____
 Address: _____
 Phone Number: _____
 Defect description: _____

Signature for acceptance of the above mentioned warranty conditions.

AVITA Corporation
 9F, No.78, Sec.1, Kwong-Fu Rd.
 Son-Chung District
 24158 New Taipei City - Taiwan

EC REP - Representative:
 Weikang Ltd., Suite B, 29 Harley Street,
 LONDON W1G 9QR, ENGLAND, United Kingdom

Distributed by:
 Medel International Srl
 Via Villagione 26 - 20156 Milano / Italy